

Open Export Utility

Many Hospitals want the endoscopy reports to be available on their Electronic Patient Record (EPR) or other similar external systems.

Unisoft's generic Open Export Utility allows reports created by Unisoft's Endoscopy Reporting Tools to be exported automatically to these external systems in a wide variety of fully configurable formats.

The export formats include Health Level Seven International (HL7), Hypertext Mark up Language (HTML), Extensible Markup Language (XML), portable document format (PDF), Comma Separated Variable (CSV), straight forward text (TXT), or a specific format for Anglia Healthcare's ICE

If you would like a demonstration or further information please call Unisoft on 0208 367 2103

